

CLEVER Semi Compliant
Tableau de Compliance
Table compliance

Pressure (bar)	1.5 mm	2.0 mm	2.25mm	2.5 mm	2.75mm	3.0mm	3.25mm	3.5mm	4.0mm	4.5mm
4	1.45	1.92	2.20	2.46	2.68	2.92	3.17	3.47	3.92	4.38
6	1.50	2.05	2.28	2.54	2.76	3.04	3.27	3.54	4.04	4.50
8	1.53	2.11	2.31	2.59	2.82	3.09	3.34	3.66	4.14	4.61
10	1.55	2.14	2.37	2.66	2.88	3.15	3.40	3.71	4.22	4.69
12	1.58	2.18	2.41	2.74	2.94	3.21	3.45	3.78	4.29	4.75
14	1.60	2.22	2.45	2.82	2.99	3.25	3.49	3.88	4.36	4.80
16	1.62	2.26	2.49	2.86	3.03	3.29	3.53	3.97	4.44	4.86
18	1.64	2.30	2.53	2.90	3.07	3.34	3.57	4.04	4.54	4.93
20	1.68	2.38	2.58	2.98	3.12	3.40	3.62	4.15	4.65	*

Gamme du produit
Product range

References	∅ Proximal/Distal Balloon (mm)	Longueur ballon (mm)
CV-150-10	1.5	10
CV-150-15	1.5	15
CV-150-20	1.5	20
CV-200-10	2.0	10
CV-200-15	2.0	15
CV-200-20	2.0	20
CV-200-25	2.0	25
CV-225-10	2.25	10
CV-225-15	2.25	15
CV-225-20	2.25	20
CV-225-25	2.25	25
CV-225-30	2.25	30
CV-225-35	2.25	35
CV-225-40	2.25	40
CV-250-10	2.5	10
CV-250-15	2.5	15
CV-250-20	2.5	20
CV-250-25	2.5	25
CV-250-30	2.5	30
CV-250-35	2.5	35

CLEVER Semi Compliant

References	∅ Proximal/Distal Ballon (mm) ∅ Proximal/Distal Balloon (mm)	Longueur ballon (mm) Balloon Length (mm)
CV-250-40	2.5	35
CV-275-10	2.75	40
CV-275-15	2.75	10
CV-275-20	2.75	15
CV-275-25	2.75	20
CV-275-30	2.75	25
CV-275-35	2.75	30
CV-275-40	2.75	35
CV-300-10	3.0	40
CV-300-15	3.0	10
CV-300-20	3.0	15
CV-300-25	3.0	20
CV-300-30	3.0	25
CV-300-35	3.0	30
CV-300-40	3.0	35
CV-325-10	3.25	40
CV-325-15	3.25	10
CV-325-20	3.25	15
CV-325-25	3.25	20
CV-325-30	3.25	25
CV-325-35	3.25	30
CV-325-40	3.25	35
CV-350-10	3.5	40
CV-350-15	3.5	10
CV-350-20	3.5	15
CV-350-25	3.5	20
CV-350-30	3.5	25
CV-350-35	3.5	30